## How to Use Alignment Tools

GA-00 for super pin design and A-600/700/800 for c clip design)

- 1. Drill hole in each side of negative mold.
- 2. Grease and insert rod through both holes.
- 3. Pour cast.
- 4. When plaster sets, remove rod and strip cast.
- 5. Pick fixture closest to size and adjust out to allow for the modification desired.
- 6. Remove pin from hinge, and contour distal section to cast, then proximal.
- 7. **For C Clips**: Apply clip and washer to fixture; pack voids with clay.

Note: For the 505-Colt, 510-Pony, and 660-Stallion (T only) Put extra washer on first, then hinge piece, the other hinge piece, then regular washer and clip.

**For Super Pins**: fasten hinges to both ends using the M2 screws provided and insert the ends into the cast



8. Pull plastic; when cool cut as you normally would. If plastic does not come off easily, cut plastic in half, at hinge center with a hot knife and take apart in two pieces.

DO NOT CUT PLASTIC HIGH OR LOW OF CENTER

**Enclose Pin and Clip** 

RIGHT





WRONG